



Protocol

DNA Isolation from Plant by nexttec™ 1^{-Step}

- nexttec™ cleanColumns -

Cat. No. 47N.010

Cat. No. 47N.050

Cat. No. 47N.250

Version 1.0

For research only

Principle

nexttec™ 1-Step is the easiest handling and fastest DNA purification system containing a single buffer system and a one-step DNA purification after lysis.

Proteins, detergents and low molecular weight compounds are retained by the nexttec™ sorbent. DNA passes through the nexttec™ cleanColumn during a short, one-step purification procedure.

The obtained DNA is suitable for all common enzymatic reactions (restriction digests, real-time PCR, PCR, genotyping etc.).

Kit contents

The kit contains all necessary reagents for lysis and subsequent DNA purification.

| Component | Art.No. 47N.010 | Art.No. 47N.050 | Art.No. 47N.250 |
|-------------------------------|----------------------------|----------------------------|----------------------------|
| Buffer F | 4,2 ml | 21 ml | 84 ml |
| Protease | 0.3 ml | 1.5 ml | 6 ml |
| Prep Solution | 6 ml | 20 ml | 100 ml |
| Solution H | 3 ml | 15 ml | 75 ml |
| SDS | 0.5 ml | 0.75 ml | 3 ml |
| DTT (1,4- Dithio-DL-threitol) | 0.1 ml | 0.25 ml | 1.5 ml |
| RNase A | 0.2 ml | 1 ml | 4 ml |
| nexttec™ cleanColumns | 10 | 50 | 250 |
| Waste collection tubes | 10 | 50 | 250 |
| DNA collection tubes | 10 | 50 | 250 |

nexttec™ service

To extend the application range to samples, which is difficult to lyse by the standard procedure, it is recommended to include optional components in the lysis buffer and optimize the lysis conditions. Please get in contact with service@nexttec.biz for detailed information.

Storage Conditions

During shipment all kit components are stable at room temperature. After arrival, store **Buffer F, Prep Solution, Solution H, DTT, Protease and RNase A** at **+2°C to +8°C**. **SDS** and **nexttec™ cleanColumns** are stored at **room temperature (+20°C to +25°C)**. If properly stored, see expiration date for the stability of the kit.

Safety Information

| | | | |
|-----------------|----|--------------------|------------------|
| Protease | Xn | R 36/37/38 R42/43; | S 23-24-26 36/37 |
| DTT | Xn | R 22-36/37/38; | S 26-36 |
| SDS | Xi | R 41; | S 25-26-39 |

Risk Phrases

| | |
|--------------------|---|
| R 22: | Harmful if swallowed. |
| R 42/43: | May cause sensitisation by inhalation and skin contact. |
| R 36/37/38: | Irritating to eyes, respiratory system and skin. |
| R 41: | Risk of serious damage to eyes. |

Safety Phrases

| | |
|-----------------|---|
| S 23: | Do not breathe gas/fumes/vapour/spray. |
| S 24: | Avoid contact with skin. |
| S 26: | In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. |
| S 36/37: | Wear suitable protective clothing and gloves. |
| S 39: | Wear eye/face protection. |

When working with chemicals, always wear a suitable lab coat, disposable gloves and protective goggles. For more information please consult the appropriate material safety data sheets (MSDS).

Before starting

- **Equilibrate nexttec™ cleanColumns**

| | |
|----|---|
| E1 | Add 350 µl Prep Solution to a nexttec™ cleanColumn . Incubate for at least 5 min at room temperature. Centrifuge at 350x g for 1 min to remove excess buffer. |
| E2 | Discard the Waste collection tube. Place the nexttec™ cleanColumn into a new DNA collection tube. Use equilibrated nexttec™ cleanColumns or store closed at +2°C to +8°C and use within one week . |

- **Preheat a thermomixer to 56°C**

Protocol

Lysis

A: Without homogenising (using leave pieces)

| | |
|----|---|
| L1 | Transfer 3 -10 pieces of leaf tissue (5x5mm) to a reaction tube. |
| L2 | Add 280 Buffer F, 20 µl Protease. Incubate with shaking (56°C, 1200 rpm, 30 min to overnight) in a thermomixer |

B: With homogenising

| | |
|----|--|
| L1 | Transfer 50 mg plant leaf tissue to a reaction tube. |
| L2 | Add 185 µl Solution H, 3 µl DTT and 12 µl RNase* and homogenise tissue in a bead mill or by grinding with a pestle. |
| L3 | Add 70 µl Buffer F, 20 µl Protease and 10 µl SDS*. Incubate with shaking (56°C, 1200 rpm, 30 min) in a thermomixer. |
| L4 | Centrifuge the lysate (10,000x g, 1 min). Use clear supernatant for DNA purification. |

C: Plant seeds

| | |
|----|---|
| L1 | Transfer up to 20 mg seed meal, grist or crude squashed seeds to a reaction tube. |
| L2 | Add 280 µl Buffer F, 20 µl Protease and 3 µl DTT*. Incubate with shaking (56°C, 1200 rpm, 30 min) in a thermomixer |
| L3 | Centrifuge the lysate at 10,000x g for 1 min Use the clear supernatant for DNA purification. |

*For Pre-Mixes see Technical Section.

Purification of DNA from lysis A,B and C

| | |
|---|---|
| P | Transfer 100 µl of the lysate to an equilibrated nexttec™ cleanColumn. Incubate for 3 min at room temperature. Centrifuge at 700x g for 1 min. |
| | The eluate contains the purified DNA!! |

Technical Section

Preparation of Lysis buffers (Pre-Mixes)

- **A: without homogenising (using leave pieces)**

| | | | | |
|--|-------------------------|----------|----------------|----------------|
| LF | Lysis Buffer LF: | 1 sample | <50 samples* | >50 samples* |
| | Buffer F | 280 µl | 280 µl x (n+3) | 280 µl x (n+5) |
| | Protease | 20 µl | 20 µl x (n+3) | 20 µl x (n+5) |
| Mix by vortexing. Add 300 µl of Buffer LF to each sample (L2). The Lysis Buffer LF is stable for 1 working day if stored at +2°C to +8°C . | | | | |

- **B: with homogenising,**

| | | | | |
|--|--------------------------|----------|----------------|----------------|
| LF1 | Lysis Buffer LF1: | 1 sample | <50 samples* | >50 samples* |
| | Solution H | 185 µl | 185 µl x (n+3) | 185 µl x (n+5) |
| | DTT | 3 µl | 3 µl x (n+3) | 3 µl x (n+5) |
| | RNase | 12 µl | 12 µl x (n+3) | 12 µl x (n+5) |
| Mix by vortexing. Add 200 µl of Buffer LF1 to each sample (L2). The Lysis Buffer LF1 is stable for 1 working day if stored at +2°C to +8°C . | | | | |
| LF2 | Lysis Buffer LF2: | 1 sample | <50 samples* | >50 samples* |
| | Buffer F | 70 µl | 70 µl x (n+3) | 70 µl x (n+5) |
| | Protease | 20 µl | 20 µl x (n+3) | 20 µl x (n+5) |
| | SDS | 10 µl | 10 µl x (n+3) | 10 µl x (n+5) |
| Mix by vortexing. Add 100 µl of Buffer LF2 to each sample (L3). The Lysis Buffer LF2 is stable for 1 working day if stored at +2°C to +8°C . | | | | |

- **C: for plant seeds**

| | | | | |
|--|-------------------------|----------|----------------|----------------|
| LF | Lysis Buffer LF: | 1 sample | <50 samples* | >50 samples* |
| | Buffer F | 280 µl | 280 µl x (n+3) | 280 µl x (n+5) |
| | Protease | 20 µl | 20 µl x (n+3) | 20 µl x (n+5) |
| | DTT | 3 µl | 3 µl x (n+3) | 3 µl x (n+5) |
| Mix by vortexing. Add 300 µl of Buffer LF1 to each sample (L2). The Lysis Buffer LF1 is stable for 1 working day if stored at +2°C to +8°C . | | | | |

*n= samples [e.g. 22 samples: Buffer F: 280 µl x (22+3)]

- **Determination of DNA concentration in nexttec™ 1-Step DNA preparations**

We recommend determining the DNA concentration:

- Using the fluorescent dye Picogreen® or similar.

- Comparing the fluorescence intensity of DNA bands of unknown concentration with standards, e.g. in ethidium bromide stained agarose gels.

Please notice:

The use of absorption measurement at 260nm (A_{260}) in a spectrophotometer (e.g. NanoDrop®) for determination of DNA concentration is system related not recommend.

For details and possible workarounds for your specific application please contact:

service@nexttec.biz.

- **Centrifugation**

For centrifugation at 350x g or 700x g use the settings for relative centrifugal force (RCF) of your centrifuge. Alternatively measure the distance of the nexttec™ cleanColumn to the centre of your rotor and calculate the necessary rotations per minute (e.g. rpm = $299.07 \times \sqrt{350/r}$; r=radius in cm)

Product Use Restriction

nexttec™ 1-Step DNA Isolation Kit components were developed, designed and sold **for research purposes only**. They are suitable for in vitro uses only. No claim or representation is intended for use to identify any specific organism or of clinical use.

It is the responsibility of the user to verify the use of the nexttec™ 1-Step DNA Isolation Kit for a specific application as the performance characteristic of this kit has not been verified to a specific organism.

Troubleshooting, FAQ and Special Applications

Product claims are subject to change. Therefore, please, visit our website or contact our technical service team for troubleshooting guide, up-to-date protocols and latest applications on nexttec™ 1-Step products.



XpressBio

Tel: 301.228.2444

Fax: (301) 560.6570