



UltraSense
***Mycoplasma* PCR Detection Kit**

Instruction Manual



Cat. No. CSHD5(sample), CDHD25, CSHD50 and CSHD100

Research Use Only. Not for Use in Diagnostic Procedures.

UltraSense PCR Detection Kit

Introduction

UltraSense Mycoplasma PCR Detection Kit utilizes the PCR, which is established as the method of choice for highest sensitivity in the detection of mycoplasma, ureaplasma and acholeplasma contamination in cell culture and other cell culture derived biologicals. The primers are specific to the highly conserved 16S ribosomal RNA sequences in the mycoplasma genome. This allows for detection of all mycoplasma species including *Acholeplasma laidlawii*, *M. arginini*, *M. fermentans*, *M. gallisepticum*, *M. genitalium*, *M. hominis*, *M. hyorhinitis*, *M. orale*, *M. pneumoniae*, *Spiroplasma citri* and *U. urealyticum*. All mycoplasma species including 8 genus and about 209 species can be detected simultaneously. The “European Pharmacopoeia” and Guideline of KFDA recommend checking for unspecific detection of *Clostridium*, *Lactobacillus*, and *Streptococcus*. The 16S ribosomal RNA of other bacteria such as *E.coli*, *Clostridium*, *Lactobacillus*, *Sterptococcus*, and plant and animal cells is not amplified.



Kit Specificity

The UltraSense Mycoplasma PCR Detection Kit detects *Mycoplasma* species simply, reliably, and rapidly. To detect the presence of these microorganisms, the assay uses the polymerase chain reaction (PCR) to amplify a target unique to a wide variety of mycoplasmas. The kit can detect more than 209 different *Mycoplasma species*, including *Acholeplasma laidlawii* and *Spiroplasma citri*. The kit does not detect other genera or cell-line DNA.

Kit Sensitivity

The sensitivity of the PCR using this kit is less than 10fg of the target DNA per reaction. Sensitivity of the assay in real culture samples depends on the quality of the sample preparation.

Preventing Template Cross-Contamination

Take precautions to minimize the potential for carryover of nucleic acids from one experiment to the next. Use separate work areas and pipettors for template preparation and PCR setup steps. Use the PCR grade water delivered with the kit, aerosol-preventive filter tips and gloves.

The 2×QPCR master mix contains dUTP instead of dTTP. When dUTP replaces dTTP in PCR amplification, UNG treatment (Uracil-N-glycosylase, not provided in this kit) can prevent the subsequent reamplification of dU-containing PCR products. UNG acts on single- and double-stranded dU-containing DNA by hydrolysis of uracil-glycosidic bonds at dU-containing DNA sites. When this strategy is put to use, carryover contamination will be eliminated while template DNA (DNA containing T) will be left intact.

Materials Provided

Materials Provided	Quantity			
	CSHD5 (Sample)	CSHD25	CSHD50	CSHD100
2x PCR Master Mix (Blue Cap)	75 $\mu\ell$	375 $\mu\ell$	750 $\mu\ell$	1.5ml
Primer Mix (Red Cap)	10 $\mu\ell$	50 $\mu\ell$	100 $\mu\ell$	200 $\mu\ell$
Positive Control DNA (Yellow Cap)	3 $\mu\ell$	13 $\mu\ell$	25 $\mu\ell$	50 $\mu\ell$
DNase Free Water (White Cap)	60 $\mu\ell$	300 $\mu\ell$	600 $\mu\ell$	1.2ml

Storage Conditions

Upon receipt, store at -20 °C. Full activity is guaranteed for 12 months.

Test Protocol

Prepare the template (Sample)

Samples should be derived from cultures which are at 90-100% confluence. Penicillin and streptomycin in the culture media do not inhibit mycoplasma or affect test sensitivity. To avoid false positive results, we recommend the use of the PCR grade water delivered with the kit, aerosol-preventive filter tips and gloves.

1. High contamination : Heat-inactivation of the sample material

The templates for the PCR analysis are prepared by direct heating of the samples (the cell culture supernatant or the biological material).

- 1) 150 μ l liquid supernatant of the sample is transferred into a sterile reaction tube.
- 2) The supernatant is centrifuged (5 minutes, 1,000 x g) to sediment cell debris.
- 3) 100ul of the supernatant is transferred into a sterile reaction tube.
- 4) The supernatant is incubated at 95°C for 10 minutes
- 5) The supernatant is used as template in the PCR. If the template contains PCR inhibition materials, the DNA can be purified with a commercial extraction kit.

2. Low contamination : Enrichment of mycoplasma by centrifugation

- 1) 1.2 ml liquid supernatant of the sample is transferred into a sterile 1.5ml reaction tube and centrifuged (5 minutes, 1,000 x g) to sediment cell debris.
- 2) 1ml of the supernatant is transferred into a sterile 1.5ml reaction tube.
- 3) The supernatant is centrifuged (10 minutes, 13,000 x g) to sediment mycoplasma particles.
- 4) The supernatant is rejected and the pellet is suspended into 50ul buffer (10 mM Tris, pH 8.0).
- 5) The suspended pellet should be vortexed and finally heated up to 95°C for 10 min.
- 6) The heated suspension can be stored at -20°C for a period of one year. Repeated freezing and defrosting, or storage in the refrigerator for longer than 12 hours should be avoided.
- 7) If the template contains PCR inhibition materials, the DNA can be purified with a commercial extraction kit.

3. Genomic DNA extraction

- 1) Prepare $5 \times 10^6 \sim 1 \times 10^7$ cells in 10ml culture
- 2) DNA was isolated using a commercial kit, QIAamp® DNA Mini (Qiagen, Valencia, CA) or equivalent products of it following the procedure provided by the vendor.
- 3) The concentration of genomic DNA was determined by UV260 measurement.
- 4) Take 5 μ L in final 50 μ L DNA solution as template for PCR reaction.

Prepare for PCR

1. Prepare the set of reactions listed in the following table. These include two types of control reactions: 1) **positive control** reaction(s) containing Mycoplasma positive control template DNA, and 2) a **negative control** (no template).

Reaction Components	Sample Reaction (Sample, μ l)	Control Reactions	
		Positive Control (Control DNA)	Negative Control (Water)
2X Master Mix	15	15	15
Primer mix	2	2	2
Template	1~10	-	-
Control DNA	-	1	-
H ₂ O (up to 30)			
Final volume	30	30	30

2. Set up the PCR instrument to run the PCR cycling (amplification) program specified below.

Steps		Temp(°C)	Time
Pre Heat		94	5 min
PCR	40 Cycles	Denature	94
		Anneal	60
		Extend	72

3. Apply 10~20 μ l each of PCR products to the gel electrophoresis.

Results

1. When mycoplasma contamination exists, a band with around 250-270bp appears. An internal DNA band with around 700bp means the right performance of PCR reaction.




Fig. 1. Lane 1, 100bp DNA ladder; lane 2, 10pg *M. pneumoniae* DNA; lane 3, 1pg *M. pneumoniae* DNA; lane 4, 100fg *M. pneumoniae* DNA; lane 5, 10fg *M. pneumoniae* DNA; lane 6, 1fg *M. pneumoniae* DNA; lane 7, 8, PCR reagent control (negative control)

Note:

- 1) Recommend to perform one negative control without sample and one positive control reaction by adding 1 μ l of mycoplasma control DNA.
- 2) If the PCR reaction is inhibited by high FBS concentration, the use of genomic DNA as a template may be helpful.
- 3) PCR inhibiting substances may accumulate in the medium of hybridoma cell. In this case, the use of diluted sample or genomic DNA as a template may be helpful.

Mycoplasma species detected by *UltraSense* Mycoplasma PCR Detection Kit

A. <i>*Acholeplasma laidlawii</i>	27. <i>Mycoplasma imitans</i>
B. <i>*Mycoplasma arginini</i>	28. <i>Mycoplasma maculosum</i>
C. <i>*Mycoplasma fermentans</i>	29. <i>Mycoplasma meleagridis</i>
D. <i>*Mycoplasma gallisepticum</i>	30. <i>Mycoplasma mycoides</i>
E. <i>*Mycoplasma hyorhinis</i>	31. <i>Mycoplasma muris</i>
F. <i>*Mycoplasma orale</i>	32. <i>Mycoplasma neurolyticum</i>
G. <i>*Mycoplasma pneumoniae</i>	33. <i>Mycoplasma opalescens</i>
H. <i>*Mycoplasma synoviae</i>	34. <i>Mycoplasma penetrans</i>
I. <i>*Spiroplasma citri</i>	35. <i>Mycoplasma pirum</i>
J. <i>Mycoplasma agalactica</i>	36. <i>Mycoplasma primum</i>
K. <i>Mycoplasma alligatoris</i>	37. <i>Mycoplasma pulmonis</i>
L. <i>Mycoplasma anatis</i>	38. <i>Mycoplasma putrefaciens</i>
M. <i>Mycoplasma arthritidis.</i>	39. <i>Mycoplasma salivarium strain</i>
N. <i>Mycoplasma bovis</i>	40. <i>Mycoplasma spermatophilum</i>
O. <i>Mycoplasma bovigenitalium</i>	41. <i>Mycoplasma sp. ovine/caprine</i>
P. <i>Mycoplasma capricolum</i>	42. <i>Mycoplasma sp. putative</i>
Q. <i>Mycoplasma cloacale</i>	43. <i>Mycoplasma suis</i>
R. <i>Mycoplasma falconis</i>	44. <i>Mycoplasma sualvi</i>
S. <i>Mycoplasma faucium</i>	45. <i>Mycoplasma timone</i>
T. <i>Mycoplasma flocculare</i>	46. <i>Ureaplasma diversum</i>
U. <i>Mycoplasma gallinarum</i>	47. <i>Ureaplasma urealyticum</i>
V. <i>Mycoplasma genitalium</i>	
W. <i>Mycoplasma hominis</i>	
X. <i>Mycoplasma pneumoniae</i>	
Y. <i>Mycoplasma hyopharyngis</i>	
Z. <i>Mycoplasma hyosynoviae</i>	

**Mycoplasma species in KFDA guidance on the Mycoplasma Test*