



One-Step RT PCR Kit

Cat. No: XV0062P-F (10 rxns)

Cat. No: XV0062P (100 rxns)

Cat. No: XV0063P (500 rxns)

Introduction

XpressBio's One-Step RT PCR Kit allows for efficient cDNA synthesis and PCR in a single tube. The kit includes a PCR master mix supplied in a 2X concentration to perform standard PCR. The Master Mix contains all the reagents (except PCR primers and template) needed for running PCR reactions. In addition, a separate RT mix that comprises a balanced mixture of both RTase and RNase Inhibitor is also provided.

RT-PCR is used to amplify double-stranded DNA from single-stranded RNA templates. In the RT step, the reverse transcriptase synthesizes single-stranded DNA molecules (cDNA) complementary to the RNA template. In the first cycle of the PCR step synthesis, Taq DNA polymerase synthesizes DNA molecules complementary to the cDNA, thus generating a double-stranded DNA template. During subsequent rounds of cycling, the DNA polymerase exponentially amplifies the double-stranded DNA template.

Features

- Higher specificity, sensitivity, and yield
- Efficient thermostable Reverse Transcriptase and RNase Inhibitor providing high cDNA yields
- Unique Hot Start Taq DNA Polymerase in a mix with high-quality dNTPs
- PCR enhancers allowing sensitive low background amplification

Application

- One-step RT-PCR
- Amplification of GC-rich and complex templates

Kit Contents

Item	100 rxns ¹	500 rxns ²
One-Step PCR Mix (2X)	2 x 1.25 ml	10 x 1.25 ml
RT Mix	2 x 125 µl	10 x 125 µl
RNase-Free Water	2 ml	5 x 2 ml

¹100 rxn of 50 µl; ²500 rxn of 50 µl

Storage

Upon receipt of the kit, immediately store the components at -20 °C. Avoid repeated freezing and thawing.

Quality Control

One step RT-PCR using eukaryotic total RNA as a template.

Protocol

The following protocol is recommended for a 50 µl reaction volume.

1. Thaw kit components, template DNA, primers, and nuclease-free water on ice. Mix each solution well.
2. Set up the following reaction mixture:

Component	Volume reaction 50 µL	Final concentration
One-Step PCR Mix (2X)	25 µL	1X
RT mix	2.5 µl	1X
Forward Primer (10 µM)	2 µL	0.4 µM ⁽¹⁾
Reverse Primer (10 µM)	2 µL	0.4 µM ⁽¹⁾
RNA template	X µL	0.1–1 µg ⁽²⁾
Nuclease-Free Water to final volume of	50 µL	

⁽¹⁾ A final primer concentration of 0.4 µM for each primer is generally optimal. However, for best results, a primer titration using 0.15–0.5 µM is recommended.

⁽²⁾ For optimal performance, use 0.1–1 µg total RNA or 10–500 ng mRNA.

3. Mix reagents completely, then transfer to a thermocycler.
4. Program the appropriate PCR cycling protocol on your PCR instrument.

Suggested thermal cycling conditions:

Step	Temperature	Time	Cycles
Reverse Transcription	45-55°C	10-20 min	1
Initial activation	95°C	2 min	1
Denaturation	95°C	10 s	35-40
Annealing *	55°C	10 s	
Extension	72°C	30-60 s/kb	
Final extension	72°C	5 min	1
Hold	4°C	∞	1

*Approximately 5°C below T_M of primers.

5. Analyze the amplification product.

Note: As with all PCR reactions, conditions may need to be optimized. You may be able to adjust your PCR conditions to optimize reaction.

PRODUCT USE LIMITATION This product is developed, designed and sold exclusively for research purposes and in vitro use only. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals.