



nexttecTM

DNA isolation systems



nexttecTM 1-Step is the easiest handling and fastest DNA purification system containing a single buffer system and a 1-Step DNA purification after lysis in 4 minutes for various samples. Available in nexttecTM cleanColumns and nexttecTM cleanPlate96 format.

1 - Step DNA Isolation System

Tissue & Cell

- Mouse Tail, Fish, Buccal Swab, Drosophila, hair, Sperm, and Typifix

Bacteria

- Gram + and Gram -

Plant

- Leaf Tissue and Seeds

Blood

- Fresh or Frozen, EDTA or Heparinized Samples

Plasmid

**Purified DNA in
only 4 minutes !**

