

Express Biotech International 4650 Wedgewood Blvd Suite 103 Frederick, MD 21703 Tel: +(Email: in Web: w

+001.301.228.2444 info@xpressbio.com www.xpressbio.com

Certificate of Analysis Helicobacter muridarum PCR Positive Control

Catalog Number: PCRPC600

Date: August 2012

Volume: 250 µl, about 250 reactions

DNA Concentration: 100 ng/µl in TE buffer. Use 0.5 to 2.0 µl per PCR reaction.

Specific Activity: At least 500 copies per µl.

Storage conditions: -20° C to -70° C

Typical PCR reaction: Total volume 50 µl, 25 µl water, 5 µl 10x PCR Buffer, 2 µl 50

mM MgSO₄, 1 μ l 10 mM dNTP mix, 1 μ l 10 μ M Forward Primer, 1 μ l 10 μ M Reverse Primer, DNA 0.5 to 1 μ l, add water to 49.5 μ l, and add last 0.5 μ l *Taq* Polymerase. Mix, spin down gently, and put tubes in thermal cycler. Alternatively, a master mix including everything except some water and DNA may be used.

Typical Thermal Cycler Program:

1.	94°C	2 min
2.	94°C	15 sec
3.	55°C	30 sec
4.	0.5°/sec	to 68°C
5.	68°C	2 min
6.	Go to 2	35 times
7.	72°C	5 min
8.	4°C	forever
9.	End	

PCR Primers:

FOR 5'-ATG ACA AAA AAA TAT TCT TTC ACA AAA CTA TTC ATT GG-3'

REV 5'-TTA TTT TAG ATT CCA TTT AAC TGC TAA ATC ATC AAT AGT C-3'

References: Comparative Medicine 54:4, August 2004, page 384. and Feng, S. et. al.

Clinical and Diagnostic Laboratory 12:4, pages 531-536, 2005.

This product is distributed for laboratory research use only. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.