

## Certificate of Analysis

### Amplified Human Liver cDNA Library

**Catalog No.:** 10106-A

**Tissue description:** Liver tissue from one 64 year old male donor

**Vector:** Express 1

**Cells used for Electroporation:** DH10BF' *E. coli* bacterial cells, T1 phage resistant

**Diluent:** 80% SOC:20% Glycerol

**Volume per Tube:** 1 ml

**Number of Tubes:** 1

**Titer:**  $5.0 \times 10^9$  cfu/ml

**Average insert size:** 1.9 kb

**Percent recombinants:** 100%

**Primer(s) for sequencing:**

M13FWD	GTA AAA CGA CGG CCA GT
M13REV	CAG GAA ACA GCT ATG AC

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.