

Express Biotech Int'l. 503 Gateway Drive West Thurmont, MD 21788 USA
 Tel:
 301.228.2444

 Fax:
 301.560.6570

 Email:
 info@xpressbio.com

 Web:
 www.xpressbio.com

Certificate of Analysis

Human Brain cDNA Library

Catalog No.: 10101-A

Tissue description: Brain tissue from three 23-27 year old male donors

Vector: pExpress 1

Cells used for Electroporation: DH10B E. coli bacterial cells, T1 phage resistant

Diluent: 80% SOC:20% Glycerol

Volume per Tube: 1 ml

Number of Tubes: 1

Titer: at least 5×10^9 cfu/ml

Average insert size: 2.1 kb

Percent recombinants: 98%

Primer(s) for sequencing: M13FWD M13REV GTA AAA CGA CGG CCA GT CAG GAA ACA GCT ATG AC

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established. **XpressBio provides innovative state-of-the-art- cDNA libraries for life science research.**