



Certificate of Analysis

U87M6 Human Glioma in pLenti6/V5-D-TOPO Standard cDNA Library

Date: December 3, 2010

Catalog No.: CSL1159, Standard cDNA Library

Tissue description: mRNA from the human tumor cell line, U87-M6.

Vector: pLenti6/V5-D-TOPO (*Eco* RI)

Cells used for Electroporation: ElectroMAX Stbl4 (grow better in LB made from components)

Diluent: 80% SOC:20% Glycerol

Number of 2059 Tubes: 5

Total volume: 8.0 ml

Titer: 8.0×10^5 cfu/ml

Average insert size: 1.6 kb

Clone (24) size range of the QC: 0.4 to 3.2 kb

Recombinants: 92%

Storage conditions: -80° C

Primer(s) for sequencing:	M13FWD	GTA AAA CGA CGG CCA GT
	M13REV	CAG GAA ACA GCT ATG AC

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.