



PRODUCT SPECIFICATION SHEET

Name:	Simian T-cell Lymphotropic Virus Type I & II (STLV-I & II) Chimeric Recombinant Antigen (HTLV Antigen is used as a Surrogate Marker)
Catalog No.:	SA832C
Lot No.:	
Specificity:	gp21 & gp46
Format:	Liquid
Purification:	98% & 95% by SDS-PAGE
Expression System:	E.Coli
Molecular Weight:	24KDa & 22.5KDa
Concentration:	1mg/ml
Supply in:	25mM Tris-C1, 0.1% SDS, pH 8.3
Storage:	Short term: 2-8°C; Long term: -20°C

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.