



Express Biotech International
P.O. Box 458
Thurmont, MD 21788-3203

Tel: +001.301.228.2444
Email: info@xpressbio.com
Web: www.xpressbio.com

CERTIFICATE OF ANALYSIS

DEPC-Treated Water

Catalog No: DTW-100A – DTW-100D

Lot No:

Format: Vials, Bottles

Product Size: 5 x 1ml, 500ml, 5 x 100ml, 1L

Purification: Autoclaved, Membrane-Filtered

Sample Treatment: DEPC-Treated

Purity or Quality Grade: RNase-Free

DNase-Free

Molecular Biology Grade

Shipping Condition: Room Temperature

Expiration: 3 years from date of shipment

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.