Express Biotech International 4650 Wedgwood Blvd Suite 103 Frederick, MD 21703 Tel: 301.228.2444 Fax: 301.560.6570 Email: info@xpressbio.com Web: www.xpressbio.com

Certificate of Analysis

Product Name: Peroxidase-conjugated Affinity Purified Goat Anti-Mouse IgG, (H+L)

Catalog #: 595-414

Lot #:

Size: 12 ml

Expiration Date:

Antibody Concentration: 1X for ELISA Procedures, prepared from lyophilized concentrated

conjugate

Peroxidase: Purified from Horseradish Root

Buffer: 0.01M Sodium Phosphate, 0.25M NaCl, pH7.6

Stabilizer: 15 mg/ml Bovine Serum Albumin (IgG-Free, Protease-Free)

Stability: 3 Years from Date of manufacture

Shipping and Storage: 2-8°C

Purity: This product was purified from antisera by immunoaffinity chromatography using antigens couples to agarose beads.

Antibody Specification: Based on immunoelectrophoresis and/or ELISA, the antibody reacts with whole molecule mouse IgG. It also reacts with the light chains of other mouse immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. The antibody has been tested by ELISA and/or solid-phase adsorbed to ensure minimal cross-reaction with human, bovine, horse, rabbit, and swine serum proteins, but it may cross-react with immunoglobulins from other species.

Preservative: None

Country of Origin: USA

This product is for in vitro RESEARCH USE ONLY, not for diagnostic or therapeutic use. This product is not a medical device.