



XpressBio
PO BOX 576
Frederick, MD 21705

Tel: 301.228.2444
Fax: 301.560.6570
Email: info@xpressbio.com
Web: www.xpressbio.com

Certificate of Analysis

NICHD_XGC_tropInt_63

Xenopus tropicalis Small Intestine
cDNA Library

Catalog number: XGC_tropInt_63

Tissue description: Total RNA isolated from small intestine tissue of wild type stage 63 *X. tropicalis* embryos.

Vector: pCS111 (Not I—Sma I)

Cells used for Electroporation: DH10B *E. coli* bacterial cells, T1 phage resistant

Diluent: 80% SOC:20%Glycerol

Number of Tubes: 1

Total Volume: 1 ml

Titer: 5×10^9 cfu/ml

Average insert size: 2.1 kb

Recombinants: 96%

Storage conditions: -80° C

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established. Express Genomics provides innovative state-of-the-art cDNA libraries for life science research.