



Express Biotech International
201 East Patrick Street, #576
Frederick, MD 21701

Tel: +001.301.228.2444
Email: info@expressbiotech.com
Web: www.expressbiotech.com

CERTIFICATE OF ANALYSIS

ANTIGEN: **LDHV**

Catalog #: AG055

Lot #:

Tissue Control Antigen: **LDHV Control Antigen**

Catalog # AG055C

Lot #

Strain: SF9-like

Volume in each vial:

Produced in: High Five insect cells

Protein Concentration: Fraction 1-4 : 0.30 mg/ml

Fraction 5-9 : 0.4 mg/ml

Buffer: Native Elution Buffer

Store frozen at -70°C or below

Expiration: 2 years from date of shipment

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been detergent treated. However, as with all biological materials, it should be handled as potentially infectious.