



Express Biotech International
503 Gateway Drive West
Thurmont, MD 21788-3203

Tel: +001.301.228.2444
Email: info@expressbiotech.com
Web: www.expressbiotech.com

CERTIFICATE OF ANALYSIS

ANTIGEN: **Theiler's Murine Encephalomyelitis Virus**

Catalog #: AG103

Lot #:

Tissue Control Antigen: **Theiler's Murine Encephalomyelitis Virus TC**

Catalog # AG103C

Lot #

Strain: Poliovirus BeAn 8386

Volume in each vial:

Produced in: BHK-21 Cells

Protein Concentration:

Buffer: 50 mM Glycine

Detergents:

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

Suggested dilution: Dilute 125 µL into 25 mL of coating buffer, then add 80-100 µL to each ELISA well

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been detergent treated. However, as with all biological materials, it should be handled as potentially infectious.