

CERTIFICATE OF ANALYSIS

Antigen: **Pneumonia Virus of Mice**

Catalog #: AG106

Lot #

Tissue Control Antigen: **Pneumonia Virus of Mice TC**

Catalog #: AG106C

Lot #

Strain: No. 15

Produced in: BHK-21 Cells

Protein Concentration:

Buffer: 50 mM Glycine

Detergents:

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been photo-chemically and detergent inactivated. However, as with all biological materials, it should be handled as potentially infectious.