



Express Biotech International  
503 Gateway Drive West  
Thurmont, MD 21788-3203

Tel: +001.301.228.2444  
Email: info@expressbiotech.com  
Web: www.expressbiotech.com

## CERTIFICATE OF ANALYSIS

ANTIGEN: **Epizootic Diarrhea of Infant Mice** Catalog #: AG112 Lot #:

Tissue Control Antigen: **Epizootic Diarrhea of Infant Mice TC** Catalog # AG112C Lot #

Strain: Rotavirus EW

Volume in each vial:

Produced in: MA-104 Cells

Protein Concentration:

Buffer: 50 mM Glycine

Detergents:

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

Suggested dilution: Dilute 125  $\mu$ L into 25 mL of coating buffer, then add 80-100  $\mu$ L to each ELISA well

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been detergent treated. However, as with all biological materials, it should be handled as potentially infectious.