



Express Biotech International
503 Gateway Drive West
Thurmont, MD 21788-3203

Tel: +001.301.228.2444
Email: info@expressbiotech.com
Web: www.expressbiotech.com

CERTIFICATE OF ANALYSIS

ANTIGEN: **Rat Coronavirus/Sialodacryoadenitis**

Catalog #: AG115

Lot #:

Tissue Control Antigen: **Rat Coronavirus/Sialodacryoadenitis TC**

Catalog # AG115C

Lot #

Strain: Coronavirus

Volume in each vial:

Produced in:

Protein Concentration:

Buffer: 0.75x PBS

Detergents:

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

Suggested dilution: Dilute 125 µL into 25 mL of coating buffer, then add 80-100 µL to each ELISA well

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been detergent treated. However, as with all biological materials, it should be handled as potentially infectious.