



Express Biotech International  
503 Gateway Drive West  
Thurmont, MD 21788-3203

Tel: +001.301.228.2444  
Email: info@expressbiotech.com  
Web: www.expressbiotech.com

## CERTIFICATE OF ANALYSIS

ANTIGEN: **Mouse Adenovirus FL**      Catalog #: AG120FL      Lot #:

Tissue Control Antigen: **Mouse Adenovirus FL TC**      Catalog # AG120FLC      Lot #

Strain: Florida

Volume in each vial:

Produced in: NIH 3T3 Cells

Protein Concentration: 1mg/ml

Buffer: 50 mM Glycine

Detergents:

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

Suggested dilution: Dilute 125  $\mu$ L into 25 mL of coating buffer, then add 80-100  $\mu$ L to each ELISA well

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been detergent treated. However, as with all biological materials, it should be handled as potentially infectious.