



Express Biotech International
503 Gateway Drive West
Thurmont, MD 21788-3203

Tel: +001.301.228.2444
Email: info@expressbiotech.com
Web: www.expressbiotech.com

CERTIFICATE OF ANALYSIS

ANTIGEN: **Encephalitozoon cuniculi** Catalog #: AG121 Lot #:

Tissue Control Antigen: **Encephalitozoon cuniculi TC** Catalog # AG121C Lot #

Strain: Sporozoan Levaditi

Volume in each vial:

Produced in: BS-C-1 Cells

Protein Concentration:

Buffer: 0.75x PBS

Detergents:

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

Suggested dilution: Dilute 125 µL into 25 mL of coating buffer, then add 80-100 µL to each ELISA well

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been detergent treated. However, as with all biological materials, it should be handled as potentially infectious.