

CERTIFICATE OF ANALYSIS

EIA ANTIGEN: **Influenza Type A Virus** Lot# A/NY090805

Catalog Number: AG123

Strain: A/NY/55/04 X PR8 X-157 R-H3

Produced in: Embryonated eggs

Hemagglutination Titration: 1:1026

Buffer: Allantoic fluid

Detergents: None

Store frozen at -70°C or below

For *in vitro* use only

This material has been heat inactivated. However, as with all biological materials, it should be handled as potentially infectious.