



Express Biotech International
503 Gateway Drive West
Thurmont, MD 21788-3203

Tel: +001.301.228.2444
Email: info@expressbiotech.com
Web: www.expressbiotech.com

CERTIFICATE OF ANALYSIS

ANTIGEN: **Parainfluenza Virus 3** Catalog #: AG138 Lot #:

Tissue Control Antigen: **Parainfluenza Virus 3 TC** Catalog # AG138C Lot #

Strain: Paramyxovirus TS-9

Volume in each vial:

Produced in: BHK-21 Cells

Protein Concentration:

Buffer: 0.75x PBS

Detergents:

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

Suggested dilution: Dilute 125 μ L into 25 mL of coating buffer, then add 80-100 μ L to each ELISA well

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been detergent treated. However, as with all biological materials, it should be handled as potentially infectious.