



Express Biotech International  
503 Gateway Drive West  
Thurmont, MD 21788-3203

Tel: +001.301.228.2444  
Email: info@expressbiotech.com  
Web: www.expressbiotech.com

## CERTIFICATE OF ANALYSIS

ANTIGEN: **Mouse Parvovirus (MPV 1 & MPV 2)** Catalog #: AG140 Lot #:

Tissue Control Antigen: **Mouse Parvovirus (MPV 1 & MPV 2) TC** Catalog # AG140C Lot #

Strain: recombinant VP2

Volume in each vial:

Produced in: E. coli

Protein Concentration: 1mg/ml

Buffer: 0.75x PBS

Detergents:

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

Suggested dilution: Dilute 125  $\mu$ L into 25 mL of coating buffer, then add 80-100  $\mu$ L to each ELISA well

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been detergent treated. However, as with all biological materials, it should be handled as potentially infectious.