



Express Biotech International
503 Gateway Drive West
Thurmont, MD 21788-3203

Tel: +001.301.228.2444
Email: info@expressbiotech.com
Web: www.expressbiotech.com

CERTIFICATE OF ANALYSIS

ANTIGEN: **Measles (Rubeola)**

Catalog No.: AG171

Lot No.:

Date of Manufacture: November 11, 2011

Strain: Edmonston

Produced in: E-6 Cells

Protein Concentration: 1.0mg/ml

Buffer: 50 mM Glycine, pH 9.0

Detergents: 0.1% CHAPS, <0.1%SDS

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been detergent treated. However, as with all biological materials, it should be handled as potentially infectious.