



Express Biotech International  
503 Gateway Drive West  
Thurmont, MD 21788-3203

Tel: +001.301.228.2444  
Email: info@expressbiotech.com  
Web: www.expressbiotech.com

## **CERTIFICATE OF ANALYSIS**

ANTIGEN: **Canine Adenovirus Type 1 Antigen** Catalog #: AG940 Lot #:

Tissue Control Antigen: **Canine Adenovirus Type 1 Antigen TC** Catalog # AG940C Lot #

Sterotype: 1

Volume in each vial:

Produced in: A-72 Cells

Protein Concentration: 1mg/ml

Buffer: 0.05M Glycine, pH 9

Detergents: <0.1% CHAPS

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

Suggested dilution: Dilute 125 µL into 25 mL of coating buffer, then add 80-100 µL to each ELISA well

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been photochemically inactivated. However, as with all biological materials, it should be handled as potentially infectious.