



Express Biotech International  
5 South Wisner Street  
Frederick, MD 21701

Tel: +001.301.228.2444  
Email: info@expressbiotech.com  
Web: www.expressbiotech.com

## CERTIFICATE OF ANALYSIS

ANTIGEN: **Feline Coronavirus (FIP) Antigen**

Catalog #: AG950 Lot #:

Tissue Control Antigen: **Feline Coronavirus (FIP) Antigen TC** Catalog # AG950C Lot #:

Strain: WSU 79-1146

Volume in each vial:

Produced in: CRFK Cells

Protein Concentration: 1mg/ml

Buffer: 0.05M Glycine, pH 9

Detergents: <0.1% CHAPS, <0.01% SDS

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

Suggested dilution: Dilute 125 µL into 25 mL of coating buffer, then add 80-100 µL to each ELISA well

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been photochemically inactivated. However, as with all biological materials, it should be handled as potentially infectious.