



Express Biotech International
503 Gateway Drive West
Thurmont, MD 21788-3203

Tel: +001.301.228.2444
Email: info@expressbiotech.com
Web: www.expressbiotech.com

CERTIFICATE OF ANALYSIS

ANTIGEN: **Feline Parvovirus (FRLV) Antigen** Catalog #: AG970 Lot #:

Tissue Control Antigen: **Feline Parvovirus (FRLV) Antigen TC** Catalog # AG970C Lot #

Strain: Cornell

Volume in each vial:

Produced in: CRFK Cells

Protein Concentration: 1mg/ml

Buffer: 50mM Glycine Buffer, pH 9, 140 mM NaCl

Detergents: <0.1% CHAPS, <0.01% SDS

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

Suggested dilution: Dilute 125 µL into 25 mL of coating buffer, then add 80-100 µL to each ELISA well

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been photochemically inactivated. However, as with all biological materials, it should be handled as potentially infectious.