



Express Biotech International
503 Gateway Drive West
Thurmont, MD 21788-3203

Tel: +001.301.228.2444
Email: info@expressbiotech.com
Web: www.expressbiotech.com

CERTIFICATE OF ANALYSIS

ANTIGEN: **Feline Herpesvirus (FVR) Antigen** Catalog #: AG980 Lot #:

Tissue Control Antigen: **Feline Herpesvirus (FVR) Antigen TC** Catalog # AG980C Lot #

Strain: C-27

Volume in each vial:

Produced in: CRFK Cells

Protein Concentration: 1mg/ml

Buffer: 0.05M Glycine, pH 9

Detergents: <0.1% CHAPS, <0.01% SDS

Store frozen at -70°C or below

Expiration: 3 years from date of shipment

Suggested dilution: Dilute 125 μ L into 25 mL of coating buffer, then add 80-100 μ L to each ELISA well

For *in vitro* use only. For Research Use Only. Not for Diagnostic or Therapeutic Purposes.

This material has been photochemically inactivated. However, as with all biological materials, it should be handled as potentially infectious.