



Express Biotech Int'l.  
503 Gateway Drive West  
Thurmont, MD 21788 USA

Tel: 301.228.2444  
Fax: 301.560.6570  
Email: info@xpressbio.com  
Web: www.xpressbio.com

## Certificate of Analysis

### NICHD\_XGC\_tropInt\_60 *Xenopus tropicalis* Small Intestine cDNA Library

**Date:** April 4, 2007

**Catalog number:** XGC\_tropInt60

**Tissue description:** Total RNA isolated from small intestine tissue of wild type stage 60 *X. tropicalis* embryos.

**Vector:** pCS111 (Not I—Sma I)

**Cells used for Electroporation:** DH10B *E. coli* bacterial cells, T1 phage resistant

**Diluent:** 80% SOC:20% Glycerol

**Number of Tubes:** 1

**Total Volume:** 1 ml

**Titer:**  $5 \times 10^9$  cfu/ml

**Average insert size:** 1.9 kb

**Recombinants:** 96%

**Storage conditions:** -80° C

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established. XpressBio provides innovative state-of-the-art cDNA libraries for life science research.