



Certificate of Analysis

Name: SRV Antigen

Catalog No.: SA830C

Lot No.:

Cell Line for Propagation: HuT 78

Buffer: 10mM Tris, 150 mM NaCl, 1 mM EDTA, pH 7.5

Preparation of Lysate: SIV purified virus adjusted to a protein concentration of approximately 1.0 mg/ml and Triton X-100 added to a final concentration of 1.0% in TNE buffer.

Protein Concentration: 1.1 mg/ml

Shipping and Storage: This product is shipped frozen on dry ice. Store at -70°C upon receipt. Avoid multiple freeze-thaw cycles as product degradation may result.

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.