



CERTIFICATE OF ANALYSIS

Name: Simian Virus 40 (SV40) Antigen

Catalog #: SA855C

Lot #:

Format: Liquid

Purification: Whole Virus/Cell lysate

Expression System: BSC-I Cells

Concentration: 1mg/ml

Supply In: 0.7 TX PRZ, 0.01% SDS, pH 7.4

Storage: Short term: 2-8°C ; Long Term: -20°C to -70°C

This material has been detergent inactivated. However, as with all biological materials, it should be handled as potentially infectious. This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.