



Express Biotech International
P.O. BOX 458
Thurmont, MD 21788 USA
Tel: 301-228-2444 Fax: 301-560-6570
Toll Free: 888-562-8914
www.xpressbio.com info@xpressbio.com

PRODUCT SPECIFICATION SHEET

Name:	Anti-Hemoglobin Monoclonal Antibody
Catalog No.:	XPFHEM-0182
Description:	DEAE Chromatography Purified
Lot No:	
Expiry Date:	
Format:	Liquid
Purity:	>95% by SDS-PAGE
Concentration:	15.22mg/ml
Supply in:	PBS (8mM Na ₂ HPO ₄ , 150mM NaCl, 2mM KH ₂ PO ₄ , 3mM KCl), PH 7.4 with 0.05% Na ₃ N
Applications:	ELISA and Lateral flow test
Storage:	Short term: 2-8°C ; Long term: -20°C Avoid multiple freeze-thaw cycles.
Unit:	Can be packed on the basis of your request

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.