

## PRODUCT SPECIFICATION SHEET

<b>Name:</b>	Anti-Hemoglobin Monoclonal Antibody
<b>Catalog No.:</b>	XPFAHEM-0183
<b>Description:</b>	DEAE Chromatography Purified
<b>Lot No:</b>	
<b>Expiry Date:</b>	
<b>Format:</b>	Liquid
<b>Purity:</b>	>95% by SDS-PAGE
<b>Concentration:</b>	13.82mg/ml
<b>Supply in:</b>	PBS (8mM Na <sub>2</sub> HPO <sub>4</sub> , 150mM NaCl, 2mM KH <sub>2</sub> PO <sub>4</sub> , 3mM KCl), PH 7.4 with 0.05% Na <sub>3</sub> N
<b>Applications:</b>	ELISA and Lateral flow test
<b>Storage:</b>	Short term: 2-8°C ; Long term: -20°C  Avoid multiple freeze-thaw cycles.
<b>Unit:</b>	Can be packed on the basis of your request

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.