



Express Biotech International

P.O.BOX 458
Thurmont, MD 21788 USA
Tel: 301-228-2444 Fax: 301-560-6570
Toll Free: 888-562-8914
www.xpressbio.com info@xpressbio.com

PRODUCT SPECIFICATION SHEET

| | |
|---------------------------|--------------------------------------------|
| Name: | Hepatitis C Virus Recombinant Antigen |
| Catalog No.: | XPFC-0014 |
| Specificity: | NS3 (1b genotype Japan) |
| Format: | Liquid |
| Expression System: | E.coli |
| Molecular Weight: | 29KDa |
| Concentration: | >0.5mg/ml |
| Supply in: | 50mM Tris, 0.1%SDS, pH8.3 |
| Applications: | ELISA &WB |
| Storage: | Short term: 2-8°C; Long term: -20°C |
| Unit: | Can be packed on the basis of your request |

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and



Express Biotech International

P.O.BOX 458

Thurmont, MD 21788 USA

Tel: 301-228-2444 Fax: 301-560-6570

Toll Free: 888-562-8914

www.xpressbio.com

info@xpressbio.com

control of the material, including its disposal, in accordance with all regulations.