



Express Biotech International
P.O.BOX 458
Thurmont, MD 21788 USA
Tel: 301-228-2444 Fax: 301-560-6570
Toll Free: 888-562-8914
www.xpressbio.com info@xpressbio.com

PRODUCT SPECIFICATION SHEET

Name: Monoclonal Antibody to Human IgG (Anti Human-IgG)
Conjugate

Catalog No.: XPFIGGP-102

Immunogen: Human IgG

Lot No.:

Format: Liquid

Purification: >90% pure.

Host Animal: Mouse, Hybridization of Sp2/0 myeloma cells with spleen cells from Balb/c mice.

Source: Ascites

Buffer: 10mMPBS, pH 7.4 and 50% glycerin

Applications: ELISA (titer need to be determined)

Recommendation: For application in indirect EIA test, this anti-human IgG Mab peroxidase conjugate is recommended for conjugating.

Storage: stored at -20°C for long term and 4°C for short term

Unit: Can be packed on the basis of your request.

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.