



Express Biotech International
P.O.BOX 458
Thurmont, MD 21788 USA
Tel: 301-228-2444 Fax: 301-560-6570
Toll Free: 888-562-8914
www.xpressbio.com info@xpressbio.com

PRODUCT SPECIFICATION SHEET

Name:	Tuberculosis Recombinant Antigen
Catalog No.:	XPFT-0040
Lot No.:	
Specificity:	16.3KD
Format:	Liquid
Purity:	95% by SDS-PAGE
Expression System:	E.coli
Molecular Weight:	16.6KDa
Concentration:	1.84mg/ml
Supply in:	50mM Tris, 0.1%SDS, pH8.3
Applications:	ELISA, rapid test &WB (titer need to be determined)
Storage:	Short term: 2-8°C ; Long term: -20°C
Unit:	Can be packed on the basis of your request.

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.