

Certificate of Analysis

Amplified Human Thymus cDNA Library

Catalog No.: 10113-A

Tissue description: Thymus tissue from one 21 year old male donor

Vector: Express 1

Cells used for Electroporation: DH10BF' *E. coli* bacterial cells, T1 phage resistant

Diluent: 80% SOC:20% Glycerol

Volume per Tube: 1 ml

Number of Tubes: 1

Titer: 5.0×10^9 cfu/ml

Average insert size: 2.0 kb

Percent recombinants: 96%

Storage conditions: -70° C

Primer(s) for sequencing:	M13FWD	GTA AAA CGA CGG CCA GT
	M13REV	CAG GAA ACA GCT ATG AC

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.