

## Certificate of Analysis

### Amplified Mouse Thyroid cDNA Library

**Catalog No.:** 10115-A

**Tissue description:** Thyroid tissue from male and female Black Swiss/Taconic mice 4-5 weeks old.

**Vector:** pExpress 1 (Not I—Eco RV)

**Cells used for Electroporation:** DH10B *E. coli* bacterial cells, T1 phage resistant

**Diluent:** 80% SOC:20%Glycerol

**Number of Tubes:** 1

**Volume per tube:** 1.0 ml

**Titer:**  $5.0 \times 10^9$  cfu/ml

**Average insert size:** 2.1 kb

**Recombinants:** 100 %

**Storage conditions:** -70° C

**Primer(s) for sequencing:**

M13FWD	GTA AAA CGA CGG CCA GT
M13REV	CAG GAA ACA GCT ATG AC

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established. Express Genomics is dedicated to providing innovative, high quality state-of-the-art cDNA libraries for life science research.