

Certificate of Analysis

Rat Spleen Primary cDNA Library

Date: May 4, 2004

Catalog No.: 20201-A

Tissue description: Total RNA purified from rat spleen.

Vector: pExpress 1 (Not I–Eco RV)

Cells used for Electroporation: DH10B *E. coli* bacterial cells, T1 phage resistant

Diluent: 80% SOC:20% Glycerol

Volume per Tube: 1 ml

Number of Tubes: 1

Titer: at least 5×10^9 cfu/ml

Average insert size: 1.7 kb

Percent Recombinants: 96%

Storage conditions: -80° C

Primer(s) for sequencing: M13FWD GTA AAA CGA CGG CCA GT
M13REV CAG GAA ACA GCT ATG AC

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.