



Express Biotech Int'l.
503 Gateway Drive West
Thurmont, MD 21788 USA

Tel: 301.228.2444
Fax: 301.560.6570
Email: info@xpressbio.com
Web: www.xpressbio.com

Certificate of Analysis

NICHD_XGC_olfb
Xenopus laevis
Embryonic Olfactory Bulb
cDNA Library

Date: April 4, 2007

Catalog No.: XGC_olfb

Tissue description: Total RNA isolated from *Xenopus laevis* embryo olfactory bulb.

Vector: pCS-111 (Not I—Sma I)

Cells used for Electroporation: DH10B *E. coli* bacterial cells, T1 phage resistant

Diluent: 80% SOC:20% Glycerol

Number of Tubes: 1

Volume per tube: 1 ml

Titer: 5×10^9 cfu/ml

Average insert size: 1.25 kb

Recombinants: 98%

Storage conditions: -80° C

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established. XpressBio provides innovative state-of-the-art cDNA libraries for life science research.