



Express Biotech Int'l.
503 Gateway Drive West
Thurmont, MD 21788 USA

Tel: 301.228.2444
Fax: 301.560.6570
Email: info@xpressbio.com
Web: www.xpressbio.com

Certificate of Analysis

NICHD_XGC_panc *Xenopus laevis* Pancreas cDNA Library

Date: September 1, 2007

Catalog No.: XGC_panc

Tissue description: Total RNA harvested from the pancreas of three adult *X. laevis*.

Vector: pCS111 (Not I—Sma I)

Cells used for Electroporation: DH10B *E. coli* bacterial cells, T1 phage resistant

Diluent: 80% SOC:20% Glycerol

Number of Tubes: 1

Volume: 1 ml

Titer: 1×10^{10} cfu/ml

Average insert size: 1.7 kb

Recombinants: 98%

Storage conditions: -80° C

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.

XpressBio provides innovative state-of-the-art cDNA libraries for life science research.