



Express Biotech Int'l.
503 Gateway Drive West
Thurmont, MD 21788 USA

Tel: 301.228.2444
Fax: 301.560.6570
Email: info@xpressbio.com
Web: www.xpressbio.com

Certificate of Analysis

NICHD_XGC_skin_m *Xenopus laevis* EmbryonicSkin cDNA Library

Date: April 4, 2007

Catalog number: XGC_skin(m)

Tissue description: Total RNA isolated from male/female wild type *Xenopus laevis* tadpoles, stage NF55, dorsal back skin tissue.

Vector: pCS111 (Not I—Sma I)

Cells used for Electroporation: DH10B *E. coli* bacterial cells, T1 phage resistant

Diluent: 80% SOC:20% Glycerol

Number of Tubes: 1

Volume per tube: 1 ml

Titer: 5×10^9 cfu/ml

Average insert size: 1.2 kb

Recombinants: 96%

Storage conditions: -80° C

This product is distributed for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established. XpressBio provides innovative state-of-the-art cDNA libraries for life science research.