

Which Chromatrap® kit is for you?

Are you working with insect samples?

Yes

Do you need NGS library reagents?

Yes

Chromatrap
Drosophila
UniqSeq
(CT-500276/CT-500277)
ProA/G

No

Chromatrap
Drosophila ChIP-Seq
(CT-500274/CT-500275)
ProA/G

No

Are you fixing your samples?

Yes

Are you using FFPE tissue?

Yes

Chromatrap FFPE
ChIP-Seq
(CT-500235/CT-500236)
ProA/G

No

Do you need NGS library reagents?

Yes

Are you using Sonication or Enzymatic digestion?

Sonication

Chromatrap
UniqSeq Sonication
(CT-500264/CT-500265)
ProA/G

Enzymatic

Chromatrap
UniqSeq Enzymatic
(CT-500266/CT-500267)
ProA/G

No

Are you using Sonication or Enzymatic digestion?

Sonication

Chromatrap ChIP-Seq
Sonication
(CT-500189/CT-500190)
ProA/G

Enzymatic

Chromatrap ChIP-Seq
Enzymatic
(CT-500191/CT-500192)
ProA/G

Chromatrap
Native
ChIP-Seq
(CT-500237/CT-500238)
ProA/G



Chromatrap

www.chromatrap.com

Distributed by:

