

XPRESS CAPTURE — COATED PLATES —

FEATURES

- Enhanced binding for sensitive assays
- Ready-to-use pre-blocked plates reduce non-specific binding
- Detergents used to lyse cells do not inhibit binding activity
- Supplied in 5-plate packs
- Available in 96 or 384-well format with black, white, or clear color options
- Stable to ship at room temperature



PRODUCTS

- Glutathione - binds fusion protein containing glutathione S-transferase (GST)
- Nickel chelate - Ni^{2+} activated surface enables metal-chelate binding of His-tagged proteins
- Streptavidin - binds to biotinylated biomolecules

