


HIV-1 REVERSE TRANSCRIPTASE ASSAY

APPLICATIONS

- Drug screening in 2 hours
- 24 hour protocol for viral conditioned media
- Qualitative measurement of HIV-1 RT with inhibitors and test articles
- Quantitative measurement of HIV-1 RT activity in samples



FEATURES

- Removable ELISA strips
- Extended shelf life
- Detects as little as 1 pg HIV-1 RT
- Other retroviruses detected include:
HIV-2, M-MuLV, AMV, SIV-1

