

HIV-1 P24 ELISA KITS


FEATURES

- Rapid Titer
- Extended shelf life
- P24 control included
- Removable ELISA Strips
- Quantitative HIV-1 P24 concentration in cell culture


SENSITIVITY

- Normal
3 to 200 pg/mL
- High
1 to 100 pg/mL
- Extended Range
100 to 4000 pg/mL

EXAMPLE STANDARD CURVE



EXAMPLE STANDARD CURVE
EXTENDED RANGE



HIV-1 PRODUCTS

- HIV-1 P24 ELISA Kit (Normal & High Sensitivity)
- HIV-1 P24 Immune Complex Dissociation Kit
- HIV-1 P24 Extended Range ELISA Kit
- HIV-1 Integrase ELISA Assay

