

AAV2 ELISA ASSAY

The AAV2 ELISA Assay is used to measure titrations for AAV2 – including AAV2 wildtype virions, assembled/intact empty AAV2 capsids, and recombinant virions.

KEY FEATURES

- **Reliable Results:** Use cell culture media instead of diluent, reduces data variation and background.
- **Long Shelf Life:** Stable positive control and 12 x 8-well breakable strips included.
- **Superior Precision:** Inter-assay CV < 10%.
- **High Sensitivity:** Minimum detection level of 2.1×10^7 capsids/ml.
- **Cost-effective:** Kit contains pre-coated 96-well ELISA plate with all reagents needed to perform the assay.
- **Rapid Results:** Total assay run time ~ 3 hours.

