

Mycoplasma species detected by BioMycoX Mycoplasma Detection Kit

1. *Acholeplasma laidlawii*
2. *Mycoplasma arginini*
3. *Mycoplasma salivarium*
4. *Mycoplasma agalactica*
5. *Mycoplasma alligatoris*
6. *Mycoplasma anatis*
7. *Mycoplasma arthritis*
8. *Mycoplasma bovis*
9. *Mycoplasma bovigenitalium*
10. *Mycoplasma capricolum*
11. *Mycoplasma cloacale*
12. *Mycoplasma falconis*
13. *Mycoplasma faucium*
14. *Mycoplasma fermentans*
15. *Mycoplasma flocculare*
16. *Mycoplasma gallinarum*
17. *Mycoplasma gallisepticum*
18. *Mycoplasma genitalium*
19. *Mycoplasma hominis*
20. *Mycoplasma hyopneumoniae*
21. *Mycoplasma hyopharyngis*
22. *Mycoplasma hyorhinis*
23. *Mycoplasma hyorhinis ATCC29052*
24. *Mycoplasma hyosynoviae*
25. *Mycoplasma imitans*
26. *Mycoplasma maculosum*
27. *Mycoplasma meleagridis*
28. *Mycoplasma mycoides*
29. *Mycoplasma muris*
30. *Mycoplasma neurolyticum*
31. *Mycoplasma opalescens*
32. *Mycoplasma orale*
33. *Mycoplasma penetrans*
34. *Mycoplasma pirum*
35. *Mycoplasma pneumonia*
36. *Mycoplasma primatum*
37. *Mycoplasma pulmonis*
38. *Mycoplasma putrefaciens*
39. *Mycoplasma salivarium*
40. *Mycoplasma salivarium strain MP1166*
41. *Mycoplasma spermatophilum*
42. *Mycoplasma sp. ovine/caprine*
43. *Mycoplasma sp. putative*
44. *Mycoplasma strain PG50*
45. *Mycoplasma suis*
46. *Mycoplasma sualvi*
47. *Mycoplasma synoviae*
48. *Mycoplasma timone*
49. *Ureaplasma diversum*
50. *Ureaplasma parvum*
51. *Ureaplasma urealyticum*