

PRODUCT SPECIFICATION SHEET

Name:	Anti-Vitamin D Mab
Catalog No.:	XPFAVD-101
Immunogen:	Vitamin D derivate-BSA
Subvariety:	IgG1
Lot No:	
Expiry Date:	
Format:	Liquid
Purity:	>90% pure
Host Animal:	Mouse
Source:	Ascites
Concentration:	4mg/ml
Supply in:	10mMPBS + 150MnNaCl, pH 7.4, 0.1% P300
Applications:	EIA
Recommended Use:	For application in ELISA
Storage:	Short term: 2-8°C ; Long term: -20°C
Unit:	Can be packed on the basis of your request

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.