



Express Biotech International
P.O.BOX 458
Thurmont, MD 21788 USA
Tel: 301-228-2444 Fax: 301-560-6570
Toll Free: 888-562-8914
www.xpressbio.com info@xpressbio.com

PRODUCT SPECIFICATION SHEET

Name:	Monoclonal Antibody to PSA HRP Conjugate
Catalog No.:	XPFPSAP-101
Immunogen:	Native PSA antigen
Lot No.:	
Format:	Liquid
Purification:	>90% pure.
Host Animal:	Mouse, Hybridization of Sp2/0 myeloma cells with spleen cells from Balb/c mice.
Source:	Ascites
Buffer:	10mMPBS, pH 7.4
Applications:	ELISA (titer need to be determined)
Recommendation:	For application in double-antibody EIA test, this anti-PSA Mab peroxidase conjugate is recommended for conjugating.
Storage:	stored at -20°C for long term and 4°C for short term
Unit:	Can be packed on the basis of your request.

This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.